










EZPass™ Suture Passers

- 1 Braided Nitinol has been designed with an oxide coating to minimize the friction transferred to soft tissue
- 2 Seven different curvatures to provide multiple options from any portal position
- 3 Sharp, low profile tips designed to ease tissue penetration and minimize tissue trauma (see back for options)
- 4 Thumb wheel allowing for quick and easy advancement of the Nitinol passing kite
- 5 Disposable instrument packaged sterile and ready for surgical use
- 6 An ergonomic handle design allows the use of a natural, relaxed, grip
- 7 Unique colors allowing for clear identification of desired style



Ordering Information

Instrumentation

Product	Description	Part Number
	Braided Nitinol Kite	110007380
	90° Up Black Handle	904050
	45° Up Gray Handle	904051
	45° Dog Leg Up Dark Blue Handle	904052
	70° Left Red Handle	904054
	70° Right Yellow Handle	904053
	30° Right Light Blue Handle	904055
	30° Left Green Handle	904056

This material is intended for health care professionals only. Distribution to any other recipient is prohibited. All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated. This material must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Check for country product clearances and reference product specific instructions for use. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and www.zimmerbiomet.com or call your local Zimmer Biomet representative.

NOT FOR DISTRIBUTION IN FRANCE

© 2018 Zimmer Biomet



0573.1-GLBL-en-REV0418 A4



Authorized Representative

Biomet UK Limited
Waterton Industrial Estate
Bridgend
CF31 3XA
United Kingdom



Legal Manufacturer

Biomet Sports Medicine
56 East Bell Drive
P.O. Box 587
Warsaw, Indiana 46581
USA

zimmerbiomet.com



CE mark on a surgical technique is not valid unless there is a CE mark on the product label.