

Shaver Blades – 130 mm

Full Radius Resector

	Design
Inner Blade	Smooth
Outer Blade	Smooth



Recommended for

stabilize meniscal rim
 full or partial excision of synovial membrane
 resection of degenerated tissue
 chondroplasty

Diameter	Item Code	Company	Product Names
4.2 mm	9181142	Smith & Nephew	-
5.5 mm	9181155	Stryker	
		Conmed Linvatec	Full Radius Resector
		Arthrex	Sabre

Shaver Blades – 130 mm

Full Radius Resector Sharp Cut

	Design
Inner Blade	Sharp
Outer Blade	Smooth



Recommended for

subacriomial decompression
notchplasty
chondroplasty

Diameter	Item Code	Company	Product Names
4.2 mm	9182742	Smith & Nephew	Bonecutter
		Stryker	Resector Cutter
		Conmed Linvatec	Ultracut
		Arthrex	Bone Cutter



Shaver Blades – 130 mm

Aggressive Full Radius

	Design
Inner Blade	Toothed
Outer Blade	Smooth



Recommended for

clear field of tissue
grab and resect frayed tissue
grab and resect cartilaginous parts and fibrillation

Diameter	Item Code	Company	Product Names
3.5 mm	9181335	Smith & Nephew	Full Radius ELITE
4.2 mm	9181342	Stryker	Aggressive Plus Cutter
4.8 mm	9181348	Conmed Linvatec	Tiger
5.5 mm	9181355	Arthrex	Dissector

Shaver Blades – 130 mm

Aggressive Full Radius Plus

	Design
Inner Blade	Toothed
Outer Blade	Toothed



Recommended for

extensive tissue debridement
resection of degenerated tissue, cartilaginous parts and fibrillation

Diameter	Item Code
4.2 mm	9183442



Company	Product Names
Smith & Nephew	Full Radius ELITE
Stryker	Aggressive Plus Cutter
Conmed Linvatec	Tiger
Arthrex	SabreTooth

Shaver Blades – 130 mm

Double Aggressive Full Radius

	Design
Inner Blade	Toothed w/ dual inner cutting windows
Outer Blade	Smooth



Recommended for

Rotator cuff prep, Meniscectomy
Chondral debridement, Synovial tissue resection

Diameter	Item Code	Company	Product Names
4.2 mm	9182342	Smith & Nephew	Synovator
5.5 mm	9182355	Stryker	Double Bite Cutter
		Conmed Linvatec	n/a
		Arthrex	Double Cut



Shaver Blades – 130 mm

Conical Double Full Radius

	Design
Inner Blade	Smooth
Outer Blade	Smooth



Recommended for

Meniscectomy

Diameter	Item Code	Company	Product Names
4.2 mm	9182442	Smith & Nephew	n/a
		Stryker	n/a
		Conmed Linvatec	n/a
		Arthrex	Torpedo

Shaver Blades – 130 mm

Aggressive Meniscus Cutter

	Design
Inner Blade	Toothed
Outer Blade	Toothed



Recommended for

Meniscectomy
Synovialectomy

Diameter	Item Code	Company	Product Names
3.5 mm	9181535	Smith & Nephew	Incisor
4.8 mm	9181548	Stryker	Tomcat Cutters
5.5 mm	9181555	Conmed Linvatec	Cuda / Great White
		Arthrex	Excalibur

Shaver Blades – 130 mm

Aggressive Meniscus Cutter Plus

	Design
Inner Blade	Toothed
Outer Blade	Toothed w/ enlarged cutting window



Recommended for

Meniscectomy
ACL stump resection

Diameter	Item Code	Company	Product Names
4.2 mm	9182142	Smith & Nephew	Incisor Plus ELITE
		Stryker	Tomcat Cutters
		Conmed Linvatec	Cuda / Great White
		Arthrex	Excalibur



Shaver Blades – 130 mm

Round Burr

Item Code	Outer Dia.	Burr Dia.	Flutes
9181042	4.2 mm	3.0 mm	8-flute
9181055	5.5 mm	4.0 mm	10-flute



Recommended for

Notchplasty
Chondroplasty

Company	Product Names
Smith & Nephew	Abrader
Stryker	Round Burr
Conmed Linvatec	Spherical Burr
Arthrex	Round Burr

Shaver Blades – 130 mm

Oval Burr

Item Code	Outer Dia.	Burr Dia.	Flutes
9181955	5.5 mm	4.0 mm	10-flute
9181965	6.5 mm	5.0 mm	10-flute



Recommended for

Subacromial decompression
Notchplasty

Company	Product Names
Smith & Nephew	Acromionizer
Stryker	Oval Burr
Conmed Linvatec	Barrel Burr
Arthrex	Oval Burr

Shaver Blades – 210 mm

Aggressive Full Radius


	Design
Inner Blade	Toothed
Outer Blade	Smooth (Cannulated)



Recommended for

clear field of tissue
grab and resect frayed tissue
grab and resect cartilaginous parts and fibrillation

Diameter	Item Code	Type
4.2 mm	9181342H	Straight
5.5 mm	9181355H	Straight
4.2 mm	9183242H	Curved
5.5 mm	9183255H	Curved



Shaver Blades – 210 mm

Aggressive Meniscus Cutter

	Design
Inner Blade	Toothed
Outer Blade	Toothed (cannulated)



Recommended for

Meniscectomy
Synovialectomy

Diameter	Item Code	Type
4.2 mm	9181342H	Straight
5.5 mm	181555H	Straight



Shaver Blades – 210 mm


Aggressive Meniscus Cutter

Item Code	Outer Dia.	Burr Dia.	Flutes
9183642H	4.2 mm	5.0 mm	10-flute



Recommended for

Femoroacetabular impingement



Shaver Blades – 70 mm

Aggressive Full Radius

	Design
Inner Blade	Toothed
Outer Blade	Smooth



Recommended for

clear field of tissue
grab and resect frayed tissue
grab and resect cartilaginous parts and fibrillation

Diameter	Item Code	Company	Product Names
2.0 mm	9181320SJ	Smith & Nephew	Full Radius ELITE
2.5 mm	9181325SJ	Stryker	Aggressive Plus Cutter
		Conmed Linvatec	Tiger
		Arthrex	Dissector

Shaver Blades – 70 mm

Round Burr

Item Code	Outer Dia.	Burr Dia.	Flutes
9181025SJ	2.5 mm	1.8 mm	6-flutes
9181035SJ	3.5 mm	2.5 mm	6-flutes



Recommended for

Small Joint Procedures

Company	Product Names
Smith & Nephew	Abrader
Stryker	Round Burr
Conmed Linvatec	Spherical Burr
Arthrex	Round Burr